

MEDICAL REGULATORS

SPECIFICALLY DESIGNED FOR MEDICAL USE

4 YEAR WARRANTY ON FUNCTION

GAS DISTINCT CASES AND GAUGE



MEDIREG - O2 - BULLNOSE PART No. 5000
PIN INDEX PART No. 5005



MEDIREG - AIR - BULLNOSE PART No. 5010
PIN INDEX PART No. 5015



DIALREG - O2 - BULLNOSE PART No. 5050
PIN INDEX PART No. 5055



FLOWREG - O2 - BULLNOSE PART No. 5150
PIN INDEX PART No. 5152

MEDIREG SPECIFICATION

- COMMON HIGH PRESSURE POD MANUFACTURED IN BRASS
- PRE SET TO APPROX. 60 PSI OUTPUT
- STANDARD MEDIREG HAS $\frac{3}{8}$ " OUTLETS
- HEAVY DUTY BULLNOSE STEM FOR ADDITIONAL SAFETY
- ALL REGULATORS SERIAL NUMBERED FOR TRACEABILITY
- GAS DISTINCT COLOUR CASES i.e. WHITE = O₂
BLACK/WHITE = AIR
- VITON 'O' RINGS ON BULLNOSE FOR INCREASED WEAR
- PIN INDEX YOKES HAVE GAS DISTINCT COLOURS
- COLOURED DISCS TO INDICATE YEAR OF MANUFACTURE

TECHNICAL SPECIFICATION

- INLET PRESSURE - 2000 PSI
- PRE SET OUTLET PRESSURE - 4 BAR (60 PSI)
- SAFETY VALVE BLOW OFF - 6 BAR (APPROX.)
- FLOW TOLERANCES - AS PER BS EN 738-1 1997
- ACCURACY +/- 10% OF INDICATED FLOW UP TO A
MAXIMUM DIFFERENCE OF 0.5 LPM


- FULL SPARES AVAILABILITY
- FULL WARRANTY ON FUNCTION FOR FOUR YEARS

THE THERAPY EQUIPMENT LTD

UNIT 1, CRANBOURNE AVENUE
CRANBOURNE INDUSTRIAL ESTATE
POTTERS BAR, HERTS. EN6 3JN

TEL: 01707 652270 FAX: 01707 652622

E-mail: sales@therapyequipment.co.uk
Website: www.therapyequipment.co.uk



0086



FM 37970