



THE THERAPY EQUIPMENT LTD



TEK 2000 EMERGENCY KITS



TEK 2000 EMERGENCY KIT - PART No. 4475

A PORTABLE SUCTION AND O₂ THERAPY KIT RUN FROM A GAS CYLINDER

SUITABLE FOR EMERGENCY SITUATIONS OR FOR USE ON THEATRE TRANSPORT TROLLEYS

KIT COMPRISES:

MINIATURE VENTURI SUCTION UNIT TO DRAW 575mm HG SUCTION FOR USE IN ON OR OFF MODE ONLY

A THERAPY EQUIPMENT DIALREG (0-15 LPM) WHICH HAS AN ADDITIONAL OUTLET TO FIX THE VENTURI SUCTION CONTROLLER

A 1/2 LITRE JAR AND CONNECTION TUBE

A CRADLE FOR THE KIT AND 1/2 LITRE JAR


(Cylinder not included)

THE THERAPY EQUIPMENT LTD

UNIT 1, CRANBOURNE AVENUE
CRANBOURNE INDUSTRIAL ESTATE
POTTERS BAR, HERTS. EN6 3JN

TEL: 01707 652270 FAX: 01707 652622

E-mail: sales@therapyequipment.co.uk
Website: www.therapyequipment.co.uk



0086



FM 37970