

THEATRE SUCTION UNIT



THEATRE SUCTION UNIT - PART No. 4760

**A STURDY 3 FOOT TROLLEY
WITH 5 CASTOR BASE
COMPLETE WITH ANGLED
HIGH SUCTION
CONTROLLER**

**THE UNIT IS POWERED
FROM THE MEDICAL
VACUUM SUPPLY AND IS
SUPPLIED WITH 5 METRES
OF VACUUM HOSE AND
INDIRECT PROBE**

**THE UNIT IS FITTED WITH
TWO 1.8 LITRE RECEIVER
JAR ASSEMBLIES AND A
CHANGE-OVER VALVE TO
SWITCH FLOWS BETWEEN
JARS EASILY**

**THE UNIT CAN BE
CUSTOMISED TO TAKE
ALTERNATIVE COLLECTION
JARS**


THEATRE SUCTION UNIT - PART No. 4760

THEATRE SUCTION UNIT - PART No. 4760

**UNIT 1, CRANBOURNE AVENUE
CRANBOURNE INDUSTRIAL ESTATE
POTTERS BAR, HERTS. EN6 3JN**

TEL: 01707 652270 FAX: 01707 652622

**E-mail: sales@therapiequipment.co.uk
Website: www.therapiequipment.co.uk**



0086



FM 37970