

REGULATORS - MULTISTAGE & ACCESSORIES



6047

Two stage adjustable Regulators ("Mickey Mouse" Type).
Adjustable between 0-150 PSI.
Predominantly for use with Air Tools.
Function Warranty: 1 Year.

MULTISTAGE REGULATORS

6040	Multi Stage Regulator - Bullnose - O2	200.00
6041	Multi Stage Regulator - Bullnose - Air	200.00
6050	Multi Stage Regulator - Pin Index - O2	260.00
6051	Multi Stage Regulator - Pin Index - Air	260.00

MULTISTAGE REGULATOR & SCHRADER VALVE

Multi Stage Regulator with Gas Distinct Schrader Valve attached.

6045	Multi Stage Regulator & Schrader Valve - Bullnose - O2	236.00
6047	Multi Stage Regulator & Schrader Valve - Bullnose - MA7	236.00
6055	Multi Stage Regulator & Schrader Valve - Pin Index - O2	296.00
6057	Multi Stage Regulator & Schrader Valve - Pin Index - MA7	296.00

NEW FOR OLD MULTISTAGE REGULATORS

We supply new equipment in exchange for your old.

5300	New for Old Multi Stage Regulator - Bullnose - O2	170.00
5301	New for Old Multi Stage Regulator - Pin Index - O2	221.00
5310	New for Old Multi Stage Regulator & Schrader - Bullnose - O2	200.60
5315	New for Old Multi Stage Regulator & Schrader - Bullnose - MA7	200.60
5311	New for Old Multi Stage Regulator & Schrader - Pin Index - O2	251.60
5316	New for Old Multi Stage Regulator & Schrader - Pin Index - MA7	251.60



6150/5
(Pack of 5)



6151




6153



6157/5
(Pack of 5)

CYLINDER KEYS

6150/5	Bullnose Cylinder Key (Spindle Key) (Pack of 5)	16.25
6151	Pin Index Cylinder Key (T Type)	17.50
6153	Ratchet Cylinder Key	22.00
6154	Combination Spanner (Open Spanner and Bullnose Key combined)	11.60
6155/5	Chain to attach key to Regulator (Pack of 5)	13.75
6157/5	Consumable Pin Index Key (Pack of 5)	13.75



6159/10 (Pack of 10)

SEALS

6159/10	Bodok Seal (Pack of 10)	9.50
6304/10	Bullnose O Ring (Pack of 10)	6.00
6099/10	Vapouriser O Ring (Pack of 10)	6.00



5701

STEPDOWN REGULATORS

Converts a 7 Bar Input Pressure to a 4 Bar Output Pressure.

5701	Stepdown Regulator Inclusive (MA7 Probe - MA4 Schrader)	123.00
------	---	--------