

SUCTION - ACCESSORIES

PIPELINE PROTECTORS (HYDROPHOBIC FILTERS)



4900

7900



4902


4903



4800



4260-05



4315/4302

Self contained Hydrophobic Filter which prevents any passage of liquid. This guarantees total protection of the hospital pipeline. Suggested change: High Usage Area - Every 3 months; Low Usage Area - Annually Or variable depending on hospital infection control policy.

4900	Bayonet Pipeline Protector (Therapy Suction) - Pack of 10	33.50
4901	Threaded Pipeline Protector (East/Old Ohmeda) - Pack of 10	62.50
4902	Small Pipeline Protector (Ohmeda Intl/Electric Suction) - Pack of 10	52.50
7900	Diamond Suction Pipeline Protector - Pack of 10	36.50
4903	Bacteriological In-Line Disc Filter - Pack of 24 (Non-Hydrophobic)	125.00

1.8 LITRE RECEIVER JARS

Manufactured in fully autoclavable Polysulphone with integral 'V' Socket. The lid is a 'pressure cooker' type bayonet fitting onto a silicone seal. Features are that it can be used as a re-usable jar or in conjunction with a Therapy Equipment Disposable Jar Liner.

4800	1.8 Litre Polysulphone Jar Assembly complete	52.50
------	--	-------

1/2 LITRE SUCTION JARS

Generally for use with Injector Suction. Does not have a float however is fully autoclavable.

4507	1/2 Litre Suction Jar	41.00
4506	1/2 Litre Suction Jar Holder (Wall Mounted)	21.00
4508	1/2 Litre Suction Jar and Wall Mounted Holder	62.00
4509	1/2 Litre Suction Jar and Rail Mounted Holder	67.00

DISPOSABLE JAR LINER

Manufactured for use with the Therapy Equipment Jar (cannot be used with other Jars) Features include:

- Non Return Valve on Patient inlet.
- Hydrophobic Filter to prevent effluent passing to controller.
- Filter Cap allowing full use of the liner before valve shut off.
- The unit is fully sealable after use.

4807	Disposable Jar Liner (Standard) - 1.8 Litre - Box of 50	330.00
------	---	--------

SUCTION ACCESSORIES

4302	Plastic Suction Tubing (30 Metre Coil)	33.00
4315	Yellow Plastic Suction Tubing (30 Metre Coil)	33.00
4260-05	Plastic V Bracket	5.50
4505	Plastic V Dovetail	7.50
4510	Abbots Bracket (Parallel)	5.50
4701-12	Microbiological Filter for Thoracic Suction Unit	5.00
4701-09	Filter Jar Only - Thoracic Suction Unit	9.00
4701-68	Complete Filter Bowl Assembly - Thoracic Suction Unit	23.50
4760-02	New/Replacement Changeover Valve for Theatre Suction Unit	52.50
4504	Suction Connector Taper	2.60
4260-60	Serres Jar Adaptor	9.50
4305	1m Yellow Tubing/Connector (Pack of 10)	35.00