



THE THERAPY EQUIPMENT LTD

# F-60B HIGH VACUUM SUCTION PUMP AC/DC

Mobile suction machine used for both Acute and Community use with a powerful double piston motor.



1

Supplied with two 2L. sterilisable jars with overflow valve. Independent use of each jar

2

ABS anti-shock case with four antistatic castors (two with brakes)

3

Complete with suction regulator and vacuum gauge

4

Supplied with disposable filter, tube sets and foot switch compatibility

## TECHNICAL DATA

Product code:	10801
Motor:	Maintenance-free Oil Free Piston Pump
Weight:	9,9 kg
Power supply:	230V/50-60Hz

Suction power:	55 l/min
Max vacuum:	- 0,80 bar - 600 mmHg
Battery duration:	> 60 min
Dimensions:	48 x 46 x 42cm

# F-60 SURGICAL SUCTION PUMP

Mobile suction unit used for both Acute and Community use with a double piston motor.



1

Supplied with two 2L. sterilisable jars with overflow valve

3

Vacuum gauge and vacuum regulator

2

ABS anti-shock case with four castors two of which with brakes

4

Supplied with disposable filter, tube sets and foot switch compatibility

## TECHNICAL DATA

Product code:	10901	Suction power:	70 l/min
Motor:	Maintenance-free Oil Free Piston Pump	Max vacuum:	- 0,80 bar - 600 mmHg
Weight:	11 kg	Dimensions:	48 x 46 x 42cm
Power supply:	230V/50-60Hz		

Therapy Equipment Ltd,  
Unit 1, Cranbourne Avenue,  
Cranbourne Industrial Estate,  
Potters Bar, Herts. EN6 3JN

Tel: 01707 652270 Fax: 01707 652622  
E-mail: sales@therapyequipment.co.uk

[www.therapyequipment.co.uk](http://www.therapyequipment.co.uk)



Distributed exclusively by  
**THERAPY EQUIPMENT LTD**

