




Therapy Equipment Ltd

DIAMOND RANGE LOW FLOW FLOWMETERS

Colour Coded for ease of Identification and Maintenance Free



GREEN BODY

0-500cc Diamond
Low Flow Flowmeter

Part No. 9583

1ST SETTING 50cc
(100cc GRADUATIONS
BETWEEN 100-500cc)



PINK BODY

0-2.5LPM Diamond
Low Flow Flowmeter

Part No. 9587

1ST SETTING 0.5LPM
(0.1LPM GRADUATIONS
BETWEEN 0.5-2.5LPM)



BEIGE BODY

0-5LPM Diamond
Low Flow Flowmeter

Part No. 9576

1ST SETTING 0.5LPM
(0.2LPM GRADUATIONS
BETWEEN 1-5LPM)

UNITS SUPPLIED WITH **7 YEAR FUNCTION AND**
ACCIDENTAL DAMAGE WARRANTY

FULLY CONFORM TO BRITISH STANDARD AND CE MARKED

Unit 1, Cranbourne Ave, Cranbourne Ind. Est., Potters Bar, Herts EN6 3JN Tel: 01707 652270



www.therapiequipment.co.uk

VERSION 3 (AUGUST 2020)



Certificate Number MD 595894