



# Therapy Equipment Ltd

Control Knob  
(over 270° Turn)  
giving control over  
full range

Red Colour coding  
denoting High Vacuum

Dual Scale  
(kPa/mmHg)  
Scale with  
'eyebrows'

Integral Metal  
Gas Distinct  
Probe

Part Number  
7701

On/Off Tap  
allowing full  
flexibility of use

Year of Manufacture  
indicator

Pipeline Protector giving  
100% guaranteed  
protection against  
pipeline contamination  
(GOES **PINK** WHEN  
ACTIVATED)

**\*\*MAINTENANCE FREE\*\***

## DIAMOND RANGE SUCTION

- CONFORMS TO BRITISH STANDARD (ISO 10079-3) & CE MARKED
- SUPPLIED WITH A LIFETIME (7 YEARS) WARRANTY AGAINST MALFUNCTION AND ACCIDENTAL DAMAGE



High Direct Diamond Suction Unit  
Part No. 7701



Rail Mounted High Diamond Suction Unit  
Part No. 7705



Low Direct Diamond Suction Unit  
Part No. 7720



Rail Mounted Low Diamond Suction Unit  
Part No. 7725



Thoracic Direct Diamond Suction Unit  
Part No. 7740



Rail Mounted Thoracic Diamond Suction Unit  
Part No. 7745

## TECHNICAL DATA

- HIGH STRENGTH PLASTIC COLOUR CODED BODY
- COLOUR CODED BACK & INDICATOR RING DENOTING TYPE OF SUCTION BEING USED:
  - RED - HIGH VACUUM (0 TO -100kPa (-750mmHg))
  - ORANGE - LOW VACUUM (0 TO -20kPa (-150mmHg))
  - GREEN - THORACIC VACUUM (0 TO -6kPa (-60cmH2O))
- GAUGE ACCURATE TO +/- 10% OF INDICATED VACUUM
- SUCTION FLOWRATE - HIGH SUCTION - IN EXCESS OF 55LPM  
LOW SUCTION - IN EXCESS OF 20LPM  
THORACIC SUCTION - IN EXCESS OF 20LPM
- SUPPLIED WITH GAS DISTINCT VACUUM PROBE, 9/16 LEFT HAND THREADED CONNECTOR, OR RAIL MOUNTED
- PIPELINE PROTECTOR FITTED AS STANDARD GIVING 100% GUARANTEED PROTECTION AGAINST PIPELINE CONTAMINATION. **FILTER GOES PINK WHEN WET**
- PRINTED RING INDICATING YEAR OF MANUFACTURE

DESCRIPTION	PART NUMBER
High Direct Diamond Suction Unit	7701
High Diamond Suction Unit with 9/16 LH Thread	7702
Rail Mounted High Diamond Suction Unit	7705
Low Direct Diamond Suction Unit	7720
Rail Mounted Low Diamond Suction Unit	7725
Thoracic Direct Diamond Suction Unit	7740
Rail Mounted Thoracic Diamond Suction Unit	7745
Diamond Range Pipeline Protector (Pack of 10)	7900

### THE THERAPY EQUIPMENT LTD

UNIT 1, CRANBOURNE AVENUE  
CRANBOURNE INDUSTRIAL ESTATE  
POTTERS BAR, HERTS. EN6 3JN

TEL: 01707 652270 FAX: 01707 652622

E-mail: sales@therapyequipment.co.uk

Website: www.therapyequipment.co.uk



Certificate Number MD 595894

VERSION 3 AUGUST 2020