



**THE THERAPY EQUIPMENT LTD**



## **TEK 2000 EMERGENCY KITS**



TEK 2000 EMERGENCY KIT - PART No. 4475

**A PORTABLE SUCTION AND O<sub>2</sub> THERAPY KIT RUN FROM A GAS CYLINDER**

**SUITABLE FOR EMERGENCY SITUATIONS OR FOR USE ON THEATRE TRANSPORT TROLLEYS**

***KIT COMPRISES:***

**MINIATURE VENTURI SUCTION UNIT TO DRAW 575mm HG SUCTION FOR USE IN ON OR OFF MODE ONLY**

**A THERAPY EQUIPMENT DIALREG (0-15 LPM) WHICH HAS AN ADDITIONAL OUTLET TO FIX THE VENTURI SUCTION CONTROLLER**

**A 1/2 LITRE JAR AND CONNECTION TUBE**

**A CRADLE FOR THE KIT AND 1/2 LITRE JAR**

**(Cylinder not included)**

**THE THERAPY EQUIPMENT LTD**

UNIT 1, CRANBOURNE AVENUE  
CRANBOURNE INDUSTRIAL ESTATE  
POTTERS BAR, HERTS. EN6 3JN

TEL: 01707 652270 FAX: 01707 652622

E-mail: [sales@therapiequipment.co.uk](mailto:sales@therapiequipment.co.uk)

Website: [www.therapiequipment.co.uk](http://www.therapiequipment.co.uk)



**VERSION 3 (AUGUST 2020)**