

THEATRE SUCTION UNIT



THEATRE SUCTION UNIT - PART No. 4760

A STURDY 3 FOOT TROLLEY WITH 5 CASTOR BASE COMPLETE WITH ANGLED HIGH SUCTION CONTROLLER

THE UNIT IS POWERED FROM THE MEDICAL VACUUM SUPPLY AND IS SUPPLIED WITH 5 METRES OF VACUUM HOSE AND INDIRECT PROBE

THE UNIT IS FITTED WITH TWO 1.8 LITRE RECEIVER JAR ASSEMBLIES AND A CHANGE-OVER VALVE TO SWITCH FLOWS BETWEEN JARS EASILY

THE UNIT CAN BE CUSTOMISED TO TAKE ALTERNATIVE COLLECTION JARS

THEATRE SUCTION UNIT

THEATRE SUCTION UNIT
UNIT 1, CRANBOURNE AVENUE
CRANBOURNE INDUSTRIAL ESTATE
POTTERS BAR, HERTS. EN6 3JN
TEL: 01707 652270 FAX: 01707 652622
E-mail: sales@therapyequipment.co.uk
Website: www.therapyequipment.co.uk



VERSION 3 (AUGUST 2020)