



THERAPY EQUIPMENT LTD

COMFORT 4 DEMAND VALVE

New improved demand valve for use with “Entonox” Gas. 4 year function and accidental damage warranty with a maximum 12 year lifespan. Fully conforms to British standards and CE Marked.



1

The Comfort Range 4 is **NOT AUTOCLAVABLE** and should be used with a single use Bacterial Filter

Suitable for ALL situations where pain relief is required (Ambulance, Emergency Services, Maternity and Birthing Pools)

2

Comfort 4 Kit with 2 Metres of PVC Hose
(Part No. 6200)


CYLINDER MODEL	Part No.
Comfort 4 Kit - Single Stage Medireg & Schrader - Pin Index - O2/N2O - Comfort 4 Demand Valve & Silicone or PVC Hose/Probe - Mouthpiece x 5 - Flo-Guard Comfort 4 Filter	
Comfort 4 Kit with 2 Metres of PVC Hose	6200
Comfort 4 Kit with 3 Metres of PVC Hose	6200/3M
Comfort 4 Kit with 4 Metres of PVC Hose	6200/4M
Comfort 4 Kit with 2 Metres of Silicone Hose	6202
Comfort 4 Kit with 3 Metres of Silicone Hose	6202/3M
Comfort 4 Kit with 4 Metres of Silicone Hose	6202/4M

PIPELINE MODEL	Part No.
Comfort 4 Kit Demand Valve - Comfort 4 Demand Valve & Silicone or PVC Hose/Probe - Mouthpiece - Flo-Guard Comfort 4 Filter	
Comfort 4 Kit Demand Valve with 2 Metres of PVC Hose	6205
Comfort 4 Kit Demand Valve with 3 Metres of PVC Hose	6205/3M
Comfort 4 Kit Demand Valve with 4 Metres of PVC Hose	6205/4M
Comfort 4 Kit Demand Valve with 2 Metres of Silicone Hose	6207
Comfort 4 Kit Demand Valve with 3 Metres of Silicone Hose	6207/3M
Comfort 4 Kit Demand Valve with 4 Metres of Silicone Hose	6207/4M

OUR PRODUCT RANGE



Medical Regulators



Flowmeters



Suction Equipment



Entonox



Portable Suction Device
and Accessories



Sundries



Therapy SOS



High/Low Air-Oxygen Mixer



www.therapyequipment.co.uk

Unit 1, Cranbourne Avenue,
Cranbourne Industrial Estate,
Potters Bar, Herts. EN6 3JN

Tel: 01707 652270 Fax: 01707 652622
E-mail: sales@therapyequipment.co.uk

